Applications of Actigraphy

By Barry T. Peterson, PhD, Sr. Mgr., Philips Respironics Clinical Affairs

About Actigraphy

Actigraphy is the measurement of motion, usually with miniature accelerometers, to monitor daily activity and sleep patterns. The devices are typically worn for several days or weeks on the wrist or the hip. The motion patterns can be displayed as actograms that show the daily activity and rest periods. The motion data can be analyzed to provide a variety of objective endpoints about the circadian patterns, the level of activity and the nighttime movements. Actigraphy has been used in many studies on various disease states.

Value of actigraphy data in studies of Depressive Disorders (depression, bipolar disorder)

Actigraphy has been used to measure the effects of various depressive disorders on sleep3, circadian patterns and daytime activity patterns^{4, 5, 7}. Although the association of circadian disorders with bipolar disorder has been well documented², more recent literature has discussed the use of sophisticated signal processing tools to show alterations in motion patterns with depression. For example, Zeitzer used functional principal component analysis of 24-hr activity patterns to differentiate among patients with Azheimer's disease who were apathetic versus those who were not⁷. Ding used similar techniques to discriminate among subjects who were healthy, mildly depressed and severely depressed¹. Krane-Gartiser used signal processing endpoints on 65 min activity samples to show that actigraphy can be used to quantify depressive moods⁵. A recent review article concludes that new analytical methods for analyzing motion data will aid in application to neurological or psychopathic disorders⁶.


Published actigraphy endpoints

Total Sleep Time, Sleep Onset Latency, Wake After Sleep Onset, Sleep Fragmentation³ Circadian endpoints (intraday variability, stability, amplitude, acrophase time)² Total activity counts every 6 hrs² Functional Principal Component Analysis^{1,4,7} Std dev. of activity, RMS successive differences, autocorrelation, entropy⁵

©2014 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

edoc SB 10/8/14 MCI 4106447



References

- Ding J, Symanzik J, Sharif A, Wang J, Duntley S, Shannon WD. Powerful actigraphy data through functional representation. Chance 2011, 24: 30-36.
- 2 Jones SH, Hare DJ, Evershek K. Actigraphic assessment of circadian activity and sleep patterns in bipolar disorder. Biopolar Disorders 2005: 7:176-186.
- 3 Kaplan, KA, Talbot LS, Gruber J, Harvey AG. Evaluating sleep in bipolar disorder: comparison between actigraphy, polysomnography, and sleep diary. Bipolar Disorders. 2012, 14: 870-879.
- 4 Kim J, Nakamura T, Kikuchi H, Sasaki T, Yamamoto Y. (2013) Co-variation of depressive mood and locomotor dynamics evaluated by ecological momentary assessment in healthy humans. PLoS ONE 8: e74979. Doi: 10.1371/journal.pone. 0074979.
- 5 Krane-Gartiser K, Henriksen TEG, Morken G, Vaaler A, Fasmer OB. Actigraphyic assessment of motor activity in acutely admitted inpatients with bipolar disorder. PLoS ONE 9(2): e89574. doi:10.1371/journal.pone.0089574.
- 6 Pan W, Song Y, Kwak S, Yoshida S, Yamamoto Y. Review article: Quantitative evaluation of the use of actigraphy for neurological and psychiatric disorders. Behavioural Neurology, Volume 2014, Article ID 897282, http://dx.doi.org/10.1155/2014/897282.
- 7 Zeitzer JM, David R, Friedman L, Mulin E, Garcia R, Wang J, Yesavage JA, Robert PH, Shannon W. Phenotyping apathy in individuals with Alzheimer disease using functional principal component analysis. Am. J. Geriatr Psychiatry 2013: 21: 391-397 (6/25).

Visit: www.actigraphy.com

Customer Service +1 541 598 3800 Toll free, US and Canada only +1 800 685 2999 Email us: minimitter@respironics@philips.com

PHILIPS

